

**TRUE**result<sup>®</sup>  
mini

BLOOD GLUCOSE MONITORING SYSTEM



## Convenience features provide competitive advantages

- No coding
- Results as fast as 4 seconds
- Tiny 0.5 microlitre sample size
- Easy-to-read display with time and date
- 500 test memory
- 7, 14 & 30 day averaging
- Strip release button
- Data management capabilities
- Proven accuracy

ISO 15197:2013  
COMPLIANT



**TRUEresult**  
mini

SMALL THINGS CAN MAKE  
A BIG DIFFERENCE

## Uses the TRUEresult Test Strips technology

- Quad electrode laser accuracy, provides precise and accurate results.
- Available in 10, 25 and 50ct test strip vials.
- Compatible with TRUEresult, TRUEresult twist, TRUEresult twist 2 and TRUEresult mini blood glucose systems.
- Exceeds the ISO 15197:2013 minimum criteria for accuracy.

All blood glucose systems are pre-set to mg/dL or mmol/L  
For more information on TRUE portfolio products, visit [www.nipro-diagnostics.eu](http://www.nipro-diagnostics.eu)